

OTHER DEVICES

EXCIMER LASERS

Company Name	Product Name	Product Description	Configuration	Working Length (cm)	Comments
Spectranetics Corporation	ELCA Laser Ablation Catheters for Coronary Intervention	0.9-mm catheter	RX	130	Spectranetics' ELCA laser ablation catheters are indicated for CTOs (crossable by a guidewire), SVGs, in-stent restenosis, PTCA dilation failures, moderately calcified lesions, ostial lesions, long lesions
		0.9- X 80-mm catheter			
		0.9-mm catheter	OTW		
		0.9- X 80-mm catheter			
		1.4-mm catheter	RX		
		1.7-mm catheter			
		1.7-mm eccentric catheter			
		2-mm catheter			
		2-mm eccentric catheter			

Recommended calibration settings: 45 Fluence, 25 Hz.

CVX-300 Excimer Laser Ablation System software upgrade required.

CTO, chronic total occlusion; OTW, over-the-wire; PTCA, percutaneous transluminal coronary angioplasty; RX, rapid exchange; SVG, saphenous vein graft.