

OTHER DEVICES

CTO DEVICES

Company Name	Product Name	Catheter Tip OD (inch)	Wire Size (inch)	Working Length (cm)	Minimum Guiding Catheter ID (F)	Comments
FlowCardia, Inc.	Crosser 14	0.043	0.014	146	6	Rapid-exchange catheter delivering high-frequency vibration for CTO recanalization
	Crosser LP	0.025	N/A		5	Low-profile catheter delivering high-frequency vibration for CTO recanalization

CTO, chronic total occlusion.