

CORONARY BALLOON-EXPANDABLE DRUG-ELUTING STENTS											
Company Name	Product Name	US FDA Indicated Use	Drug Delivered	Material Used		Polymer Type	Maximum Guidewire Size (inch)	Introducer Size (F)	Stent Diameter (mm)	Stent Length (mm)	Delivery System Length (cm)
Abbott Vascular	Xience V Everolimus-Eluting Coronary Stent System	Coronary	Everolimus	Cobalt chromium		Fluorinated copolymer	≤ 0.014	N/A	2.5, 2.75, 3, 3.5, 4.0	8, 12, 15, 18, 23, 28	143
Boston Scientific Corporation	Taxus Liberté Stent	Coronary	Paclitaxel	316L stainless steel		Translute (Styrene-b isobutylene-b styrene)	≤ 0.014	5 F compatible	2.5, 2.75, 3, 3.5, 4	8, 12, 16, 20, 24, 28, 32	OTW, 135; MR, 139
	Taxus Liberté Atom Stent								2.25		
	Taxus Liberté Long Stent								2.75, 3, 3.5, 4	38	OTW, 135; MR, 140
	Promus Stent System		Everolimus	Cobalt chromium		Fluorinated copolymer		≥ 5	2.5, 2.75, 3, 3.5, 4	8, 12, 15, 18, 23, 28	OTW, RX 143
Cordis Corporation	Cypher Sirolimus-Eluting Coronary Stent System	Coronary	Sirolimus	316L stainless steel		Nonerodible	0.014	6	2.5, 3, 3.5	8, 13, 18, 23, 28, 33	OTW
								5–6	2.25, 2.5, 2.75, 3, 3.5		RX
Medtronic, Inc.	Endeavor Sprint Zotarolimus-Eluting Coronary Stent System (RX and OTW)	Coronary	Zotarolimus	Cobalt alloy		Phosphorylcholine	≤ 0.014	≥ 5	2.5, 3, 3.5	8, 9, 12, 14, 15, 18, 24, 30	140 ± 5
MR, monorail; OTW, over the wire; RX, rapid exchange.											