

SHEATHS/WIRES/CATHETERS

STANDARD SHEATHS

Company Name	Product Name	Size (F)	Length (cm)	Guidewire Diameter (inch)	Curved or Straight	Radiopaque Tip (Yes/No)
Boston Scientific Corporation	Super Sheath	4–9	11, 25	0.035, 0.038	Straight	No
Cook Medical	Performer: Small Valve	4–8	5, 5.5, 7, 13, 25, 30, 35, 45	0.035, 0.038	Straight	With or without
Cordis Corporation	Avanti+ Catheter Sheath Introducer	4	N/A	0.035	N/A	N/A
		5–11		0.038		
	Avanti+ Transradial Catheter Sheath Introducer	4–7	11	0.021	Straight	
		4–6	23			
	Valveless/Portless Sheath Introducers	6–8	N/A	N/A	N/A	
	Long Sheath Introducers and Introducer Sets	7–8	43, 96		Straight, curved	
	Monitoring Sheaths	5–8	11, 23		N/A	
Galt Medical Corporation	Micro-Access Elite HV hemostasis valve introducer	4–7	5.5, 12	0.018	Straight	With or without
	ReDial 4 cm	5–7	4	0.018, 0.038		
Medtronic, Inc.	Input TS	5, 6, 7, 8, 9, 10, 11	7, 11, 23	N/A	Straight	No
	Input PS	5, 6, 7, 8, 9	11, 23			
Merit Medical	Prelude Sheath Introducers	4, 5, 6, 7, 8	11, 23	0.035, 0.038	Straight	With or without
	Prelude Sheath Introducers (ACT sheaths)	4.5, 5.5, 6.5, 7.5, 8.5	11			No
	Prelude Sheath Introducers (radial, etc.) (also includes 21-gauge needle with 0.018-inch wire kits and 20-gauge needle with 0.025-inch wire kits)	4–6	7, 11, 23	0.018, 0.025		
Oscor Inc.	Adelante (Peel Away)	6, 7, 7.5, 8, 8.5, 9, 9.5, 10, 10.5, 11, 12	16, 23	0.035, 0.038	Straight	No
	Adelante Sigma (Hemostic Valve)	4, 5, 6, 7, 8, 9, 10, 11	10, 23			
	Adelante Sigma Sol (Hemostic Valve)					Yes
St. Jude Medical	Ultimum	4–12, 14, 16, 18, 20, 22	5, 8.5, 12, 23, 25, 30, 35	0.035, 0.038	Straight	Yes (EV line)
	Engage	4, 5, 6, 7, 8, 9	12, 25			No

SHEATHS/WIRES/CATHETERS

STANDARD SHEATHS (CONTINUED)

Company Name	Product Name	Size (F)	Length (cm)	Guidewire Diameter (inch)	Curved or Straight	Radiopaque Tip (Yes/No)
Terumo Interventional Systems	Pinnacle Introducer Sheath for Coronary Applications	4, 5, 6, 7, 8, 9, 10, 11	10, 25 (25 does not come with mini guidewire or in 4 F)	Comes with or without 0.038 guidewire; guidewire for the 4 F is 0.035	Straight	No
	Pinnacle TIF Tip Introducer Sheath for Coronary Applications	5, 6, 7, 8	10 with 0.038-inch mini guidewire	N/A	N/A	N/A
	Pinnacle R/O II Radiopaque Marker Introducer Sheath	4, 5, 6, 7, 8, 9	6, 10, 25	Mini guidewire available for 6 and 10 cm; maximum guidewire is 0.038		Yes
	Pinnacle R/O II HiFlo	6, 7, 8	4	Comes with or without 0.035 guidewire		
	Glidesheath (indicated for radial access)	4, 5, 6	10, 16, 25	Comes with 0.021 guidewire	Spring/straight	No