# BIORESORBABLE TERMINOLOGY

# **DURABLE POLYMER**

# **PROPERTIES**

- Metal stent
- Drug embedded in permanent polymer layer
- Drug elution time varies depending on platform

# REPRESENTATIVE PRODUCTS

■ Ion (Boston Scientific Corporation),
Promus (Boston Scientific Corporation),
Resolute (Medtronic), Xience (Abbott Vascular)

# **BIORESORBABLE POLYMER**

## **PROPERTIES**

- Metal stent
- Polymer dissolves over 2–4 months
- Drug elutes as polymer dissolves

# REPRESENTATIVE PRODUCTS

■ Synergy (Boston Scientific Corporation)

# **POLYMER-FREE**

### **PROPERTIES**

- Metal stent
- Drug applied directly to stent
- Drug elutes over 1–3 months

# REPRESENTATIVE PRODUCTS

■ BioFreedom (Biosensors International)\*

\*Not available for sale in the United States.

# BIORESORBABLE SCAFFOLD (BRS, BVS)

### **PROPERTIES**

- Polymer stent
- Drug coated
- Entire stent dissolves over 2–4 years

# REPRESENTATIVE PRODUCTS

■ Absorb (Abbott Vascular)\*

\*Not available for sale in the United States.