

EMBOLIZATION COILS

Company Name	Product Name	Material Used	Coil Diameter (mm)	Coil Length (cm)	Recommended Delivery Catheter Sizes (ID)(inch)	Indicated Use
Boston Scientific Corporation ^a	VortX-18 Diamond Shape Fibered Platinum Coils	Platinum	Small secondary/large secondary: 2/3, 2/4, 2/5, 2/6	2.3, 4.1, 5.8, 8	0.021	0.018-inch fibered platinum coils are intended for arterial and venous embolizations in the peripheral vasculature
	VortX-18 Fibered Platinum Coils			2.2, 4.2, 6, 8.5		
	Straight 0.018-inch Fibered Platinum Coils		Straight	0.2, 0.5		
	Figure-8 0.018-inch Fibered Platinum Coils		2	1		
	Multi-Loop 0.018-inch Fibered Platinum Coils		4	3		
	Complex Helical 0.018-inch Fibered Platinum Coils		4, 6, 7, 8, 10, 11	2, 3, 4, 5, 6, 7		
	Multi-Loop 0.035-inch Fibered Platinum Embolization Coils		3, 4, 5, 6, 7, 9	2, 3, 4, 5, 6	0.038	For arterial and venous embolizations in the peripheral vasculature
	VortX-35 Vascular Occlusion Coils		Apex diameter/ base diameter: 2/4, 3/5, 3/6, 3/7	3, 3.5, 5.3, 6.7		
	Interlock Fibered IDC Occlusion System	Platinum-tungsten alloy	Two-dimensional helical: 3, 4, 5, 6, 8, 10, 12, 14	6, 8, 10, 12, 15, 20, 30	0.021	To obstruct or reduce rate of blood flow in the peripheral vasculature; device is not intended for neurovascular use
			VortX Diamond: small secondary/ large secondary 2/3, 2/4, 2/5, 2/6	2.3, 4.1, 5.8, 8		

^aPlease reference appropriate product Instructions for Use from Boston Scientific for a list of indications, warnings, precautions, and potential adverse effects.

OTHER DEVICES

EMBOLIZATION COILS (CONTINUED)

Company Name	Product Name	Material Used	Coil Diameter (mm)	Coil Length (cm)	Recommended Delivery Catheter Sizes (ID) (inch)	Indicated Use
Cook Medical	Nester Embolization Coils	Platinum	4, 6, 8, 10, 12	14	0.035	Arterial and venous embolization
	MicroNester Embolization Coils		3, 4, 6, 8, 10		0.018	
	Tornado Embolization Coils		5-3, 6-3, 7-3, 8-4, 8-5, 10-5	4.1–12.5	0.035	Embolization of selective vessel supply to AVMs, trauma, bleeding, and other vascular lesions
	Tornado Embolization Microcoils		3-2, 4-2, 5-2, 6-2, 7-3, 8-4, 10-4	2–14.2	0.018	
	Hilal Embolization Coils		Straight to 5	0.5–3		Embolization of selective vessel supply to AVMs and other vascular lesions of the brain, spinal cord, and spine
	MR eye Embolization Coils	Inconel	3, 5, 8, 10, 12, 14, 15, 20, 25	2, 3, 4, 5, 8, 10, 12, 15	0.035 or 0.038	Arterial and venous embolization
	Stainless Steel Embolization Coils	Stainless steel	6, 8, 10, 15, 20	8, 10, 15	0.052	
	MR eye Flipper Detachable Embolization Coils		3, 5, 6.5, 8	3, 4, 5, 6, 8, 10, 12	0.041	Arterial and venous vessel embolization in the peripheral vasculature
Terumo Interventional Systems	Azur Peripheral HydroCoil Pushable .018-inch System	Platinum and hydrogel	2, 3, 4, 5, 6, 8, 10	2, 4, 6, 10, 14	0.022	Arterial and venous vessel embolization in the peripheral vasculature
	Azur Peripheral HydroCoil Pushable .035-inch System		4, 5, 6, 8, 10, 15	4, 6, 10, 14	0.035	
	Azur Peripheral HydroCoil Detachable .018-inch System		3, 4, 6, 8, 10, 12, 15, 20	5, 10, 15, 20, 30	0.022–0.025	
	Azur Peripheral Framing Coil Detachable .018-inch System	Platinum	8, 10, 14, 20	20, 26, 34, 50		

AVM, arteriovenous malformation; ID, inner diameter.