

PTA BALLOONS

Company Name	Product Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Shaft Length (cm)	Sheath Sizing Size (F)	Balloon Diameters (mm)	Balloon Lengths (cm)				
Abbott Vascular	Viatrac 14 Plus Peripheral Dilatation Catheter	8	14	0.014	80, 135	4	4, 4.5	1.5, 2, 3, 4				
						5	5, 5.5, 6, 6.5, 7					
	Agiltrac 0.035 Peripheral Dilatation Catheter			0.035	55, 80, 135	4	2, 4, 6, 8, 10					
						5						
						6	2, 4, 6, 8					
						6	10					
						5	7	2				
						6	4, 6, 8, 10					
						8	2, 4, 6, 8					
						9	2, 4, 6					
						10	2, 4					
						7	6					
	12					2, 4						
	14											
	Agiltrac 0.018 Peripheral Dilatation Catheter					8	14	0.018	80, 135	5	4, 5, 6, 7	2, 3, 4, 6
										8	2, 3, 4	
							12			6	6	
										9	2, 3, 4, 6	
		10	2, 3, 4									
		7	6									
	6	8	3, 4									
	Fox Plus	5–10	Up to 18	0.035	50, 80, 135	5–7	3, 4, 5, 6, 7, 8, 9, 10, 12	2, 3, 4, 6, 8, 10, 12				
	Fox sv	6–12	Up to 22	0.018	90, 150	4	2, 2.5, 3, 4, 5, 6	1.5, 2, 3, 4, 6, 8, 10, 12				
	FoxCross	5–10	Up to 18	0.035	50, 80, 135	5–7	3, 4, 5, 6, 7, 8, 9, 10, 12, 14	2, 3, 4, 6, 8, 10, 12				
AngioDynamics, Inc.	WorkHorse II	10	14, 16, 18	0.035	50, 75, 100	5–8 mm through a 6 F, 9–10 mm through a 7 F	5, 6, 7, 8, 9, 10	2, 4, 6				
	Profiler	7	12, 14, 16		40, 75, 120	4–8 X 4 mm through a 5 F, 8 X 6–10 mm through a 6 F	4, 5, 6, 7, 8, 9, 10	2, 4, 6, 8				

BALLOONS

PTA BALLOONS (CONTINUED)

Company Name	Product Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Shaft Length (cm)	Sheath Sizing Size (F)	Balloon Diameters (mm)	Balloon Lengths (cm)
AngioScore, Inc.	AngioSculpt Scoring Balloon Catheter	8	20 (2–2.5 mm), 18 (3 mm), 16 (3.5 mm)	0.014	137	5	2, 2.5, 3, 3.5	1, 2
			16 (4 mm), 14 (5, 6 mm)	0.018	50 (6 mm only), 90, 137	6	4, 5, 6	2, 4
Bard Peripheral Vascular, Inc.	Conquest	8	Up to 30	0.035	50, 75	6–8	5, 6, 7, 8, 9, 10, 12	2, 3, 4, 6, 8
	Atlas	4, 6	Up to 18		75, 120	7–12	12, 14, 16, 18, 20, 22, 24, 26	2, 4, 6
	Rival	6, 8	Up to 14		80, 135	5–7	4, 5, 6, 7, 8, 9, 10	2, 4, 6, 8, 10, 15
	Optiplast XT/XL	5, 6, 8	Up to 15		50, 75, 100, 120	5–8	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10
	Dorado	8	Up to 24		40, 80, 120, 135	5–7	3, 4, 5, 6, 7, 8, 9, 10	1.5, 2, 3, 4, 6, 8, 10, 12, 15, 17, 20
	Ultraverse		Up to 15	0.018	100, 120	4, 5	2, 3, 4, 5, 6, 7	2, 4, 6, 10
	VascuTrak2	8 (2–3.5 mm), 6 (4–7 mm)	12	0.014 (2–3.5), 0.018 (4–7 mm)	140	5 (2–4 mm), 6 (5 mm), 7 (6–7 mm)	2, 2.5, 3, 3.5, 4, 5, 6, 7	2, 4, 6, 8, 10, 12, 15, 20, 25, 30
Biotronik, Inc.	Pheron	7	15	0.035	80, 130	5–6	3, 4, 5, 6, 7, 8, 9, 10	2, 4, 6, 8, 10
B. Braun	Mini Ghost	Up to 10	Up to 20	0.018	40, 80, 120, 150	4	2–6	1, 2, 3, 4, 6, 8, 10
	Z-Med	Up to 6	Up to 15	0.035	80, 100, 120	5–16	2–40	2, 3, 4, 5, 6
	Impact		Up to 10		75, 120	8–14	12–25	2, 4
	CoEfficient		Up to 13	0.018	75, 135	4–7	4–12	2, 3, 4

BALLOONS

PTA BALLOONS (CONTINUED)

Company Name	Product Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Shaft Length (cm)	Sheath Sizing Size (F)	Balloon Diameters (mm)	Balloon Lengths (cm)
Boston Scientific Corporation	Symmetry Balloon Catheter with Glidex Hydrophilic Coating	6	15	0.018	90, 135	4 (1.5–4 mm X 4 cm), 5 (4 mm X 10 cm–6 mm X 4 cm)	1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6	2, 4, 10
	Symmetry Stiff Shaft with Glidex Hydrophilic Coating Balloon Catheter				90, 135, 150			
	XXL Balloon Catheter	Not provided	8 (12, 14 mm), 5 (16, 18 mm)	0.035	75, 120	7 (12, 14 mm), 8 (16, 18 mm)	12, 14, 16, 18	2, 4, 6
	Blue Max Balloon Catheter	6	20 (4–8 mm X 4 cm), 17 (8 mm X 8 cm–10 mm)		40, 75, 120	6 (4–6 mm), 7 (7–8 mm), 8/9 (9–10 mm)	4, 5, 6, 7, 8, 9, 10	2, 3, 4, 8, 10
	Ultra-thin Diamond Balloon Catheter		15 (3–8 mm), 12 (9–10 mm), 10 (12 mm)		40, 75, 120, 135, 150	5 (3–6 mm), 6 (7–8 mm), 7 (9–12 mm)	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10
	Ultra-thin SDS Balloon Catheter		12		50, 75, 90, 135	5 (up to 6 X 4), 6 (up to 8 X 4), 7 (9–10 mm)	4, 5, 6, 7, 8, 9, 10	1.5, 2, 3, 4, 6, 8
	Synergy Balloon Catheter	6 (3–10 mm), 4 (12 mm)	12 (3–8 mm), 10 (9–10 mm), 6 (12 mm)		50, 75, 90, 135, 150	5 (up to 6 X 4), 6 (up to 8 X 4), 7 (9–12 mm)	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10
	Sterling Monorail PTA Balloon Dilatation Catheter	6	14 (3–7 mm), 12 (8 mm)	0.018, 0.014	80, 135	4 (up to 8 X 4), 5 (up to 8 X 6)	3, 3.5, 4, 4.5, 5, 5.5, 6, 6.5, 7, 8	1, 1.5, 2, 3, 4, 6
	Sterling Over-the-Wire Balloon Dilatation Catheter		14 (4–7 mm), 12 (8 mm), 10 (9–10 mm)		40, 80, 135	4 (up to 8 X 4), 5 (up to 10 X 4), 6 (up to 10 X 8)	4, 5, 6, 7, 8, 9, 10	2, 3, 4, 6, 8, 10
	Sterling ES PTA Balloon Dilatation Catheter		14 (1.5–3 mm), 12 (3.5–4 mm)	0.014	142–146	4	1.5, 2, 2.5, 3, 3.5, 4	2, 3, 4

BALLOONS

PTA BALLOONS (CONTINUED)

Company Name	Product Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Shaft Length (cm)	Sheath Sizing Size (F)	Balloon Diameters (mm)	Balloon Lengths (cm)
Cook Medical	Advance 14LP Low Profile PTA Balloon Catheter	8	12	0.014	170	4	2, 2.5, 3, 4	2, 4, 6, 8, 12, 16, 20
	Advance 18LP Low Profile PTA Balloon Catheter		14	0.018	80, 135	4 (up to 7 X 10), 5 (8 mm)	3, 4, 5, 6, 7, 8	2, 4, 6, 8, 10
	Advance 35LP Low Profile PTA Balloon Catheter	5–10	Up to 15	0.035		5 (up to 8 X 8), 6 (up to 10 X 8), 7 (12 x 14)	3, 4, 5, 6, 7, 8, 9, 10, 12	2, 3, 4, 6, 8, 10, 12, 14
	ATB Advance	5	15		40, 80, 120	5–7	4, 5, 6, 7, 8, 9, 10, 12, 14	2, 3, 4, 5, 6, 7, 8
Cordis Corporation	Savvy Long PTA Dilatation Catheter	6	15, 13, 12	0.018	120, 150	4–5	2, 2.5, 3, 3.5, 4, 5, 6	12, 15, 22
	Sleek RX PTA Dilatation Catheter		15, 16	0.014	150, 155	4	2, 2.5, 3, 3.5, 4	4, 8, 10, 12, 15, 22
	Opta Pro PTA Dilatation Catheter	6, 4 (12 mm)	10, 6 (12 mm)	0.035	80, 110, 135	5–7	3–10, 12	1, 1.5, 2–4, 6, 8, 10
	Powerflex P3 PTA Dilatation Catheter	10, 6 (12 mm)	15, 8 (12 mm)		40, 65, 80, 110, 135	5–8	4–10, 12	1, 1.5, 2–4, 6, 8, 10
	Powerflex Extreme PTA Dilatation Catheter	15, 12 (10 mm)	20, 17 (9–10 mm)		40, 80, 120	5–9	4–10	2, 4, 6
	Maxi LD PTA Dilatation Catheter	4, 3.5 (22, 25 mm)	6, 5 (22, 25 mm)		80, 110	8–11	Iliac: 14, 15	2, 4, 6, 8
	Slalom PTA Dilatation Catheter	10 (3–6 mm), 8 (6–7 mm), 6 (8mm)	14 (3–5 mm), 12 (6–7 mm), 10 (8 mm)	0.018	80, 120, 135	4–6	3–8	2, 4
	Savvy PTA Dilatation Catheter	8 (2–4 mm), 6 (4.5–6 mm)	10		80, 120, 150	4–5	2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6	2–4, 6, 10
	Aviator Plus PTA Dilatation Catheter	10	14		0.014	142	4 (4–6 mm), 5 (7 mm)	4, 4.5, 5, 5.5, 6, 7

PTA BALLOONS (CONTINUED)

Company Name	Product Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Shaft Length (cm)	Sheath Sizing Size (F)	Balloon Diameters (mm)	Balloon Lengths (cm)
ev3 Inc.	EverCross 0.035 PTA Balloon	10, 7	20, 18, 14, 12, 11, 10	0.035	40, 80, 135	5, 6, 7	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10, 12, 15, 20
	NanoCross 0.014 PTA Balloon	10, 8, 7	14	0.014	90, 150	4	1.5, 2, 2.5, 3, 3.5, 4	2, 4, 8, 12, 15, 21
Invatec, Inc.	Amphirion Deep OTW PTA Balloon Catheter	7	14	0.014	120, 150, 210 (tapered)	4	1.5, 2, 2.5, 3, 3.5, 4	2, 4, 8, 12
	Admiral Xtreme OTW PTA Balloon Catheter	8, 6 (10, 12 mm)	Up to 14	0.035	80, 130	5–7	3, 4, 5, 6, 7, 8, 9, 10, 12	2, 4, 6, 8, 12
	Submarine Plus OTW PTA Balloon Catheter	7		0.018	90, 130	4–5	2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 7	
Medtronic, Inc.	Sprinter	8, 6 (1.5 mm)	14, 12 (1.5 mm)	N/A	N/A	138	1.5, 2, 2.25, 2.5, 2.75, 3, 3.25, 3.5, 3.75, 4	6, 10, 12, 15, 20, 25, 30
	NC Stormer	10	18				2, 2.25, 2.5, 2.75, 3, 3.25, 3.5, 3.75, 4	6, 11, 16, 21, 26, 31
	Optimus PTA Balloon Catheter	8	12 (5–8 mm), 10 (9, 10 mm)	0.035	80, 135	5 (5–8 mm), 6 (9, 10 mm)	5, 6, 7, 8, 9, 10	2, 3, 4, 6, 8, 10, 12, 15