

SHEATHS/CATHETERS

STANDARD SHEATHS

Company Name	Product Name	Size (F)	Length (cm)	Guidewire Diameter (inch)	Curved/ Straight	Radiopaque Tip (Yes/No)
Arrow International, a division of Teleflex Medical	Cath Lab/Interventional Sheaths (adult)	4–9	11	0.035	Straight	Yes
	Cath Lab/Interventional Sheaths (pediatric)	4–6	5, 7.5	0.021		
	Arrow-Flex PSI (with Arrow Raulerson Syringe and Cath-Gard Catheter Contamination Guards with and without ArrowGard)	6, 7, 8, 8.5, 9	7, 10	0.025–0.035		
B. Braun	Intradyn Arterial Introducer Sheath and Elite HV Hemostasis Introducers	3–14	7, 11, 12, 23, 24	0.018, 0.021, 0.035, 0.038	Straight	With or without
Boston Scientific Corporation	Super Sheath Introducers	4, 5, 6, 7, 8, 9	11, 25	With wire, 0.035; without wire, 0.038 compatible	Straight	No
	Super Sheath XL Introducers	10,11,12,14				
Cook Medical	Performer: Small Valve	4–8	5, 5.5, 7, 13, 25, 30, 35, 45	0.035, 0.038	Straight	With or without
	Performer: Large Valve	9–18	13, 30, 35, 40, 70, 80	0.038		
	Performer: Radial	4–6	5.5, 7, 13, 23	0.018		No
Cordis Corporation	Brite Tip Interventional Sheath	4–11	5.5, 11, 23, 35, 45, 55, 90	0.035	Straight	Yes
	Avanti Sheath Introducer		5.5, 11, 23			No
ev3 Inc.	Introducer Sheaths	6	5.5	0.035	Straight	Yes
Galt Medical Corporation	Micro-Access Elite HV Hemostasis Valve Introducer	4–7	5.5, 12	0.018	Straight	With or without
	ReDial 4 cm	5–7	4	0.018, 0.038		
Medtronic, Inc.	Input TS	5, 6, 7, 8, 9, 10, 11	7, 11, 23	0.038	Straight	Yes
	Input PS	5, 6, 7, 8, 9	11, 23			
Merit Medical	Prelude Sheath Introducers	4, 5, 6, 7, 8	11, 23	0.018, 0.025, 0.035, 0.038	Straight	With or without
	Prelude ACT Sheath Introducers	4.5, 5.5, 6.5, 7.5, 8.5	11	0.035, 0.038		Without
St. Jude Medical, Inc.	Ultimum	4–12, 14, 16, 18, 20, 22	5, 8.5, 12, 23, 25, 30, 35	0.035, 0.038	Straight	Yes, EV line
Terumo Interventional Systems	Pinnacle Introducer Sheath for Peripheral Applications	4, 5, 6, 7, 8, 9, 10, 11	10, 25 (25 does not come with mini guidewire or in 4 F)	Comes with or without 0.035 guidewire	Straight	No
	Pinnacle Introducer Sheath for Coronary Applications			Comes with or without 0.038 guidewire; the guidewire for the 4 F is 0.035		
	Pinnacle TIF Tip Introducer Sheath for Peripheral Applications	5, 6, 7, 8	10 with 0.035-inch mini guidewire	Mini guidewire included		
	Pinnacle TIF Tip Introducer Sheath for Coronary Applications		10 with 0.038-inch mini guidewire			

SHEATHS/CATHETERS

STANDARD SHEATHS (CONTINUED)

Company Name	Product Name	Size (F)	Length (cm)	Guidewire Diameter (inch)	Curved/ Straight	Radiopaque Tip (Yes/No)
Terumo Interventional Systems (continued)	Pinnacle R/O II Radiopaque Marker Introducer Sheath	4, 5, 6, 7, 8, 9	6, 10, 25	Mini guidewire included with 6 and 10 cm; maximum guidewire is 0.038	Straight	Yes
	Pinnacle R/O II HiFlo	6, 7, 8	4	Comes with or without 0.035 guidewire		
	Glidesheath (indicated for radial access)	4, 5, 6	10, 16, 25	Comes with 0.021 guidewire; 6-F size also available with 0.035 guidewire	Spring/ straight	No