

OTHER DEVICES

VASCULAR PLUGS

Company Name	Product Name	Material Used	Sizes (mm)	Delivery Cable Length (cm)	Recommended Guide Catheter Sizes (ID) (inch)	Indicated Use
AGA Medical Corporation	Amplatzer Vascular Plug	Nitinol	4–8	135	0.056	Arterial and venous embolizations in the peripheral vasculature
			10–12		0.067	
			14–16		0.088	
	Amplatzer Vascular Plug II		3–8		0.056	
			10–12		0.07	
			14–16		0.086	
			18–22		0.098	

ID, inner diameter.