

OTHER DEVICES

SNARES

Company Name	Product Name	Catheter Size/ Length (F/cm)	Working Diameters (mm)	Snare Length (cm)	Snare Type
Abbott Vascular/ Vascular Solutions, Inc.	Expro Elite	0.035 inch	5, 10, 15, 25, 35	150	Helical loop
Cook Medical	Curry Intravascular Retriever	8/100, 6.3/100	0.021 inch, 0.018 inch	150	Folds at midpoint and forms a loop when passed through the catheter; supplied sterile in peel-open pack- ages; intended for one-time use
	Dotter Intravascular Retriever	8/95	30-mm basket	7-cm basket	Helical loop basket is loaded into the hub of the catheter and expands upon emergence from the distal tip; supplied sterile in peel-open packages; intended for one-time use
	Indy OTW Vascular Retriever	8/55, 8/100	40-mm basket	Not provided	Four-loop nitinol wire snare; over-the- wire capabilities; compatible with 0.035-inch guidewire
ev3 Inc.	Amplatz Goose Neck Snare Kit	4/102	5, 10	120	Single 90° loop configuration
		16/162		65	
		6/102	15, 20, 25, 30, 35	120	
		16/223		65	
	Amplatz Goose Neck Microsnare Kit	2.3–3/150	2, 4, 7	175	
2.3–3/175		200			
Alligator Retrieval Device	Not provided	Not provided	Not provided	Not provided	
Merit Medical	EN Snare Mini Snare System	3.2/150	2–4, 4–8	175	Three interlaced nitinol loops that self-size to the vessel diameter
Vascular Solutions, Inc.	Micro Elite Snare	0.014 inch	2, 4, 7	180	Helical loop
	Expro Elite Snare	0.035 inch	5, 10, 15, 25, 35	150	