

OTHER DEVICES

SPECIALTY BALLOONS

Company Name	Product Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing (F)	Shaft Length (cm)	Balloon Diameters (mm)	Balloon Lengths (cm)
Atrium Medical Corporation	ClearWay OTW Therapeutic Perfusion Catheter	2–4	N/A	N/A	Introducer size: 7, 8	N/A	40, 80	4, 5, 6, 7, 8	1
	ClearWay Rx Therapeutic Infusion Catheter				Introducer size: 4, 5	Guide catheter compatibility: 6, 7	134	1, 1.5, 2, 2.5, 3, 3.5, 4	2, 5
Boston Scientific Corporation	PolarCath Peripheral Dilatation System	8	N/A	0.035	6, 7, 8	6 (4–5 mm), 7 (6–7 mm), 8 (8 mm)	80, 120	4–8	2, 4, 6, 8, 10
				0.014	5, 5.5, 6	5, 5.5, 6	135, 150	2, 2.5, 3, 3.5, 4, 5, 6	2, 4, 6, 10
	Peripheral Cutting Balloon (OTW)	6	10	0.018	6, 7	6 (5, 6 mm), 7 (7, 8 mm)	50, 90, 135	5–8	2
	Peripheral Cutting Balloon (Monorail)		12	0.014	6	6	140	2, 2.5, 3, 3.5, 4	1.5