

OTHER DEVICES

ATHERECTOMY DEVICES

Company Name	Product Name	Minimum Vessel Diameter (mm)	Sheath Compatibility (F)	Crossing Profile (inch)	Working Length (cm)	Tip Length (cm)	
Cardiovascular Systems Inc.	Diamondback 360° 1.25-mm Classic Crown	1.5	6	0.07 (1.8 mm)	145	1	
	Diamondback 360° 1.5-mm Classic Crown	1.75					
	Diamondback 360° 1.75-mm Classic Crown	2					
	Diamondback 360° 2-mm Classic Crown	2.25		0.079 (2 mm)		1.5	
	Diamondback 360° 2.25-mm Classic Crown	2.5	7	0.089 (2.25 mm)		2	
	Diamondback 360° 1.5-mm Solid Crown (30 micron)	1.75	6	0.070 (1.8 mm)		1	
	Diamondback 360° 1.75-mm Solid Crown (30 micron)	2					
	Diamondback 360° 2-mm Solid Crown (30 micron)	2.25		0.079 (2 mm)			1.5
	Diamondback 360° 2.25-mm Solid Crown (30 micron)	2.5	7	0.089 (2.25 mm)		2	
ev3 Inc.	TurboHawk LS-C	3.5–7	8	0.105 (2.7 mm)	110	6	
	TurboHawk LX-C				113	9	
	TurboHawk SL-M		7, 8		110	6	
	TurboHawk LX-M				113	9	
	SilverHawk LS-M	4.5–7	7, 8	0.105 (2.7 mm)	110	6	
	SilverHawk LX-M	4.5–6.5			113	9	
	SilverHawk MS-M	3.5–5	110		6		
	SilverHawk SX	3–3.5	7	0.095 (2.4 mm)	136	4.3	
	SilverHawk SXL					7.2	
	SilverHawk SS+			0.09 (2.3 mm)	135	2.6	
	SilverHawk EXL	2.3	6			0.08 (2 mm)	6
	SilverHawk ES+	2–2.5				0.075 (1.9 mm)	2.2
	SilverHawk DS	1.5–2				0.073 (1.9 mm), 0.062 (1.57 mm) without wire	2.6
	Pathway Medical Technologies, Inc.	Jetstream G3	3	7		0.083 (2.1 mm)	135