

Indications/Contraindications:

Intended use: The Pounce™ Thrombectomy System is intended for the non-surgical removal of thrombi and emboli from the peripheral arterial vasculature.

Contraindications:

- The device is not intended for venous applications
- The device is not intended to be deployed in native vessels smaller than 3.5 mm or larger than 6 mm
- The device is not intended for peripheral vasculature dilatation
- The device is not for coronary or neurovascular use
- The device is contraindicated for use in patients who cannot receive recommended intravenous anticoagulant therapy
- The safety and effectiveness of the device has not been established in pediatric patients (<18 years of age)
- The device is not intended to be deployed in vessels with previously implanted devices

References (from pages 14-16)

1. Zeller T, Müller C, Frank U, Bürgelin K, Horn B, Roskamm H. The Straub-Rotarex thrombectomy system: initial experiences. *Rofo*. 2001;173:626-631. doi: 10.1055/s-2001-15843
2. de Athayde Soares R, Matielo MF, Neto FCB, et al. Analysis of the results of endovascular and open surgical treatment of acute limb ischemia. *J Vasc Surg*. 2019;69:843-849. doi: 10.1016/j.jvs.2018.07.056
3. Ansel GM. Presented at the 44th Charing Cross Symposium; April 22, 2022; London, United Kingdom.